


(logo)	TECHNICAL DATA SHEET	
		DATE: dd/mm/yy
		3-17AY-2017

<b>PRODUCT</b>	PERMLASTIC REGULAR PKG	
<b>REFERENCE</b>	60101	
<b>COMPANY/Manufacturer</b>	Kerr Italia	
<b>CLASSIFICATION</b>	CLASS I MEDICAL DEVICE Directive 93/42/CEE	
<b>PRODUCT DESCRIPTION</b>	Permlastic is available in three viscosities. Light Body is excellent for inlay and fixed bridge impressions. Regular Body is recommended for partial- or full-denture impressions. The denser viscosity of Heavy Body results in optimum compression of the Light Body in a double-mix or copper-band impression.	
<b>SPECIFICATIONS</b>		
<b>PACKAGING</b>	Permlastic Regular, 130 ml tube	
<b>SHELF LIFE</b>	720 days	
<b>INTENDED OF USE</b>	Elastic Impression Material	
<b>SECURITY</b>	Only profesional use	
<b>TRANSPORT AND STORAGE CONDITIONS</b>		

Sign/stamp

**KerrHawe, S.A.**  
Of. Representación en España  
Cartagena, 187-191, Entlo. Of. 8  
08013 Barcelona  
Tel. 93 245 29 82